




#4-40 LOCKING SETSCREW (INCLUDED)
0.050in[1.3mm] HEX
2 PLACES, BOTH SIDES

SM1 CARRIAGE LINEAR DISPLACEMENT
0.079in[2.0mm]



FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION 		YiiZeeLABS	
DATE	04/16/2025		
For reference only — not intended for manufacturing.			
	MATERIAL #		REV
COPYRIGHT  YiiZeeLABS			A
VALUE IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS		APPROX WEIGHT	kg