





Cross section A-A



4-M4*0.7

Tighten with a 2mm wrench



model	L0
CSJ30-RMS	RMS (0.8"-36)
CSJ30-C	1.00"-32UNF (C-Mount)

DRAWING PROJECTION 		YiiZeeLABS	
DATE	04/16/2025		
For reference only — not intended for manufacturing.			
	MATERIAL #		REV
COPYRIGHT  YiiZeeLABS			A
VALUE IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS		APPROX WEIGHT	kg