



ACCEPTS $\varnothing 25.4\text{mm}$ [$\varnothing 1.00\text{in}$] OPTICS
 $\varnothing 24.1\text{mm}$ [$\varnothing 0.95\text{in}$] THRU TO CENTER




3/16in-100
ACTUATOR PROVIDES MIRROR SURFACE TRANSLATION
2.0mm [$5/64\text{in}$] HEX

3/16-100 ACTUATOR WITH REMOVABLE KNOB
2.0mm [$5/64\text{in}$] HEX
2 PLACES


NYLON TIP SET SCREW
FOR SECURING OPTIC
2.0mm [$5/64\text{in}$] HEX






M6 x 1.0 TAPPED HOLE
8.9mm [0.35in] DEEP




METRIC ID MARK
M4 x 0.7 TAPPED HOLE
6.4mm [0.25in] DEEP



$\varnothing 26.3$ [$\varnothing 1.03$]mm [$\varnothing 1.04\text{in}$] THRU TO OPTIC FACE
SM1 ($\varnothing 1.035-40$) SERIES INTERNAL THREAD
0.15in [3.8mm] DEEP
2 PLACES @ 90° SEPARATION



DRAWING PROJECTION 				
DATE	2017/8/1 9			
KCD1L-M1 笼式直角调整架SM1(4-40螺纹) 装配图				
www.Oeabt.com		MATERIAL #		REV
COPYRIGHT © 2019 BY Oeabt				A
VALUE IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUNDOFF ERRORS		APPROX WEIGHT		kg