






Technical Requirements:

1. Material Specification: Stainless steel
2. Appearance Requirements: Surface must be smooth and flat, free of dents, scratches, and sharp edges. All sharp edges must be rounded. Unless otherwise specified, chamfers are C0.3.
3. Unspecified Tolerances: Conform to ANSI B4.1-1967, Class Fine (f)
4. Color: Natural (original metal finish)

DRAWING PROJECTION 		EZ Lab	
DATE	14/11/17		
OHD-3030			
		MATERIAL #	REV
COPYRIGHT  EZ Lab			A
VALUE IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			kg
		APPROX WEIGHT	