






- 技术要求:
1. 材料规格: 不锈钢
 2. 外观要求: 表面光洁, 平整, 无碰划伤锐边磨圆, 未注倒角为C0.3
 3. 未注公差; 按GB1804-f。
 4. 原色

DRAWING PROJECTION 		  Optics Experiment and Basic Teaching		
DATE	20/11/17			
OHD-4060				
www.Oeabt.com		MATERIAL #		REV
COPYRIGHT  2017 BY Oeabt				A
VALUE IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			APPROX WEIGHT	kg