




截面 A-A



| 型号 | L0 |
|---------|-------|
| PCA3-A | 30mm |
| PCA4-A | 40mm |
| PCA5-A | 50mm |
| PCA6-A | 60mm |
| PCA8-A | 80mm |
| PCA10-A | 100mm |



- 技术要求:
1. 材料规格: SUS304
 2. 外观要求: 表面光洁, 平整, 无碰划伤锐边磨圆, 未注圆角为R0.5
 3. 未注公差; 按GB1804-f。
 4. 原色

| | | | |
|--|------------|---------------------|-----|
| DRAWING PROJECTION | | <h1>YiiZeeLABS</h1> | |
| DATE | 04/16/2025 | | |
| For reference only — not intended for manufacturing. | | | |
| | | MATERIAL # | REV |
| COPYRIGHT © YiiZeeLABS | | | A |
| VALUSE IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS | | APPROX WEIGHT | kg |